

PANCREATIC ELASTASE, STOOL - Details

 About This Test

Component Results

Component	Your Value	Standard Range
Patient Value Pancreatic Elas Unit: ug E/g stool	322.1	
Interpretation Pancreatic ELAS	Normal -----ADDITIONAL INFORMATION----- ----- Reference Values For Pancreatic Elastase in Stool >200 ug Elastase/g stool = Normal 100 to 200 ug Elastase/g stool = Moderate to slight exocrine pancreatic insufficiency <100 ug Elastase/g stool = Severe exocrine pancreatic insufficiency Test Performed by: Joli Diagnostic, Inc. 2451 Wehrle Drive Williamsville, NY 14221	

General Information

Collected: 11/29/2017 1:10 PM

Resulted: 12/28/2017 1:34 PM

Ordered By: Amit Rastogi, MD

Result Status: Final result

This test result has been released by an automatic process.

MyChart® licensed from Epic Systems Corporation, © 1999 - 2015.